

# ***USER MANUAL***

## ***Vacuum Mixer***

### ***LZ-ZKJ-1***



## Description

The Vacuum mixer LZ-ZKJ-1 is used for mixing and vibrating investment materials and molded materials. As the mixing process is done in a vacuum state, the workpiece is free of bubbles, and firmness and compactness of the mold has also been improved.

## Specification

|               |                 |
|---------------|-----------------|
| Voltage       | AC220V/110V±10% |
| Powder        | 150W            |
| Bowl Capacity | 500ml           |
| Gross Weight  | 16kg            |
| Dimensions    | 32*30*58cm      |

## Method of Installation

The LZ-ZKJ-1 Vacuum Mixer is a table type machine; it will work normally combined with the air compressor. The machine should be operated following the instruction.

## Instruction

1. Unpack the package, then check if all the accessories are inside.
2. Plug in the machine.
3. Connect the machine to an air compressor.
4. Stir the material inside the vacuum bowl to pasty state.
5. Screw up the air outlet bolt. Switch on the machine and preset the **Mixing** time. Press the **Vacuum Pump** button. Connect the

vacuum bowl closely with the machine and then press ***Mixing***.

6. After finishing the mixing, switch off the machine.

### **Cautions**

1. Read all the operation steps before use to avoid damage.
2. The air pressure of the air compressor must be above 0.6mpa.
3. Do not disassemble the machine without professional instruction.
4. Please check the grounding of the machine regularly. Make sure there is a safe connection between the case of machine and the ground.
5. Make sure the investment is **less than one** half of the vacuum bowl in case there is blocking of the air circuit.

Contact Us

E-mail Address:

Dentalsalemail.com@gmail.com

Telephone Number: +86-023-68438829

Official Site:<https://www.dentalsalemail.com>