

Automatic Programmable Vacuum Porcelain Furnace

Digital-control ceramic firing equipment for dental laboratories

1190°C

MAX. TEMP.

60

PROGRAMS

<35 mmHg

FINAL VACUUM

26 kg

NET WEIGHT

PRODUCT HIGHLIGHTS

- **Programmable firing**
60 stored firing programs
- **Bilingual controls**
Chinese / English panel selection
- **Temperature adjustment**
Keypad calibration and correction
- **Automatic functions**
Auto purge and night mode
- **Motorized furnace head**
Automatic lifting and lowering
- **Cooling functions**
Cooling and multi-window display

TECHNICAL SPECIFICATIONS

Input voltage	110V / 220V (variants listed)
Rated power	1200W + 350W
Maximum temperature	1190°C
Heating rate	30-140 (unit not specified)
Holding time	0-40 minutes
Ultimate vacuum	Less than 35 mmHg
Program capacity	60
Protection rating	IPX1
Firing chamber	Diameter 8 x 4 cm
Net weight	26 kg
Shipping weight	35 kg
Package dimensions	37 x 54 x 46 cm

PRODUCT INFORMATION

Heating-rate units were not specified on the seller listing.
Confirm voltage version, pump inclusion and accessories before purchase.