



Ultrasonic Scaler VRN's New "Sonic" Era Scaler

DA-10



Built-in High-Performance Chip

Equipped with an STMicroelectronics chip, featuring advanced algorithms that automatically detect resonance frequency, preventing instantaneous impact and ensuring a comfortable scaling experience.

DA-10

Smart Touch System

Responsive waterproof touch panel with one-touch switching between **G** (scaling), **P** (periodontal), and **E** (endodontic) modes. Offers **12 power** levels for precise treatment customization.

High-Quality Handpiece

New handpiece assembled technology with a detachable design, able to use high-temperature and high-pressure sterilization.



Intelligent Feedback

Automatically adjusts power based on calculus hardness, increasing power as needed for enhanced cleaning efficiency.

SUPRAGINGIVAL SCALING



Before



After



Before



After

Supragingival Scaling

G1



Supragingival scaling, used to remove large supragingival calculus in all directions.

G2



Flat edge is used to remove heavy supragingival calculus.

G4



Used to remove all supragingival calculus and bacterial plaque and supragingival interdental areas.

Periodontal Treatment

P1



Used for subgingival calculus removal and periodontal irrigation.

