


TUOJIAN TJ2688 A1

Complete Dental Chair Unit

An integrated dental treatment unit with a DC-motor chair, instrument delivery system, water supply and suction accessories.

PRODUCT LAYOUT



Visual guide only - appearance and installed options may vary.

PRODUCT FEATURES

● DC Motor Chair

Motor-driven chair positioning

● Adjustable Headrest

Headrest position adjustment

● Instrument Delivery

Handpiece tubing and holders

● Dental Operation Light

Dual-mode light assembly

● Built-in Water System

Water bottle and water supply

● Assistant Equipment

Suction and saliva ejector

TJ2688 A1

Technical Specifications

Configuration information for the TJ2688 A1 listing

Model	TJ2688 A1
Input voltage	AC 110-240 V, 50/60 Hz (listing)
Air pressure	0.55 MPa
Water pressure	0.2-0.4 MPa
Chair height range	440-860 mm
Spittoon	Ceramic (specification field)
Instrument tray	Bottom-mounted standard; top-mounted optional
Color choices	7 color options listed
Net / gross weight	200 kg / 250 kg
Packing dimensions	136 x 95 x 112 cm

LISTED EQUIPMENT

- Electric valve control
- DC motor
- Thermostatic water supply
- Handpiece tubing, 3 sets
- Three-way syringes, 2 sets
- Saliva ejector with switch
- Suction apparatus
- Operation light, dual-mode
- X-film viewer
- Built-in water bottle
- Balance arm assembly
- Adjustable headrest
- Doctor stool, 1 set
- Foot pedal, 1 set

ORDERING / INSTALLATION NOTES

Confirm the electrical version, instrument-tray orientation, spittoon material and upholstery color before shipment. The seller listing describes both glass and ceramic spittoons; ceramic is shown above because it appears in the specification table. Installation should follow the supplied manufacturer instructions.

Contact Us

E-mail Address:

Dentalsalemall.com@gmail.com

Telephone Number: +86-023-68438829

Official Site: <https://www.dentalsalemall.com>