MODEL AX-PT1

Pressure Polymerizer

Dental Laboratory Press & Polymerization Unit

TECHNICAL DATASHEET

PRODUCT OVERVIEW

The AX-PT1 is a combined dental laboratory pressing and polymerization unit for denture processing applications. Its upright press frame and integrated control base are designed for a compact benchtop setup.

KEY FEATURES

- Combined dental press and polymerization functions
- Upright metal frame with T-shaped handle and circular pressing plates
- Mesh guard around the central work area
- Front control dial, indicator light and power switch

TECHNICAL SPECIFICATIONS

Brand / Model	AiXin AX-PT1
Power supply	AC 110 V / 220 V, 50/60 Hz
Power consumption	600 W
Motor power	6 W
Listed pressing load	2,000-3,500 kg*
Unit dimensions (H x W x D)	40 x 28 x 14 cm
Shipping dimensions (H x W x D)	46 x 34 x 25 cm
Unit weight	Approx. 10 kg

PACKAGE CONTENTS

1 x AX-PT1 polymerizer / press unit

1 x User manual

PRODUCT INFORMATION NOTE

* The published specifications use kg for a claimed pressing load (described as "rated pressure"). This is not a pressure value in bar or MPa. Confirm the nameplate, voltage version and applicable specifications with the supplier. This datasheet is not an operating manual.

dentalsale11.com